

MATERIAL SAFETY DATA SHEET according to 2001/58/EC

1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY

1.1 Identification of the substance/ preparation ENAPOST
1.2 Application of the substance / the preparation Post for aesthetic restoration
1.3 Identification of the company MICERIUM SpA Via G. Marconi, 83 16030 AVEGNO (GE) ITALIA

1.4 Emergency Info: +39 0185 7887870

2. COMPOSITION / INFORMATION ON INGREDIENTS

Chemical Characteristics: Glass fibers, 1-4 Butandiol dimethacrylate, Urethandimethacrylate, Bis-GMA

Hazardous Components: CAS-Nr Name % Reclassification
n.a.

This product doesn't contain hazardous material (for Phrases R see sez. 16)

3. HAZARDOUS IDENTIFICATION

Most important hazards n.a.
Eyes: n.a.
Skin: n.a.
Inhalation: n.a. For finishing, use respiratory protection and dust suction
Ingestion: n.a.
Exposure effects n.a.
Chronic effects: n.a.

4. FIRST-AID MEASURES

After eye contact: n.a.
After skin contact: n.a.
After inhalation: n.a.
After swallowing: Consult a doctor

5. FIRE-FIGHTING MEASURES

Extinguishing media that should not be used for safety reasons n.a.
Exposure hazard to substance / preparation n.a.
Fire or explosion hazards (conditions to avoid) n.a.
Protective equipment for fire fighters n.a.

6. ACCIDENTAL RELEASE MEASURES

Person-related safety precautions n.a.
Measures for environmental protection n.a.
Measures for cleaning/collecting: n.a.

7. HANDLING AND STORAGE

7.1 Handling:
Avoid touching posts with bare hands as body oils can compromise bond strengths.
Accidental release measures n.a.
Local and general ventilation: n.a.
Environmental precautions: n.a.
7.2 Storage
Material to be avoided n.a.
Prevention from electricity n.a.
7.3 Special uses n.a.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Additional information for technical plants n.a.

8.1 Limits of values for exposure
Components with limits of values to be supervised at the workplace:
CAS-NR Name % Art Value Unit

The product contains no relevant quantities of components with limits of values to be supervised at the workplace.

8.2 Exposure controls Necessary only during finishing
8.2.1 Professional Exposure Control
8.2.1.2 Protection of hands: Protective gloves
8.2.1.3 Protection of eyes: Safety glasses
8.2.1.4 Protection of skin: Use respiratory protection and dust suction for finishing,
8.2.2 Environmental Exposure control n.a.

=====

9. PHYSICAL AND CHEMICAL PROPERTIES

=====

9.1 General information

Form: conical rod
Color: Enamel plus UD1 shade
Odor: none

9.2 Important informations on human health, safety and environment

Information on change in the physical state	Value	Unit	Method
PH: n.a.			Ripartition coefficient: n.a.
Boiling point: n.a.			Water: n.a.
Flash point: n.a.			Content of solvents: n.a.
Danger of explosion: n.a.			Organic solvents: n.a.
Density: n.a.			
Vapor density : n.a.			Content of solids: n.a.
Viscosity: n.a.			9.3.1 Other informations
Combustion property: n.a.			Miscibility in: insoluble
Vapor pressure: n.a.			Conductivity: n.a.
Solubility in: insoluble			Melting point: n.a.
Hydrosolubility: n.a.			Gas group: n.a.
Liposolubility (solvent – fat): n.a.			Autoignition temperature: n.a.

=====

10. STABILITY AND REACTIVITY

=====

10.1 Conditions to be avoided n.a.
10.2 Materials to be avoided n.a.
10.3 Dangerous decomposition products: n.a.

=====

11. TOXICOLOGICAL INFORMATION

=====

Toxicological effects from exposure to substance and preparation n.a.
Primary irritant effect:
- on the skin/ on the eye n.a.
- inhalation n.a.
- ingestion n.a.
Sensitization of product: n.a.

=====

12. ECOLOGICAL INFORMATION

=====

General information n.a.

=====

13. DISPOSAL CONSIDERATION

=====

Dispose in accordance with Federal, State and Local regulations

=====

14. TRANSPORT INFORMATION

=====

Maritime transport (IMDG): n.a.
Land transport (ADR/RID): n.a.
Air transport (ICAO / IATA): n.a.
UN number / exact shipment denomination / packaging n.a.
Additional information

=====

15. REGULATORY INFORMATION

=====

Informations on health, safety, environment n.a.

=====

16. OTHER INFORMATION

=====

Symbols and Indication of Danger (Phrases S and R): n.a.

FOR PROFESSIONAL USE ONLY

Further informations: n.a.

Revision Summary: No revisions have been made.

The information and recommendations set forth herein (hereinafter "information") are presented in good faith and are believed to be correct as of the date hereof. Micerium Spa, however makes no representations as to the completeness of security of this information and supplies it on the condition that the persons receiving same will make their own determination as to its suitability for their purposes prior to use. In no event will Micerium Spa, be responsible for damages of any nature whatsoever resulting from the use of or reliance upon information.

NO REPRESENTATIONS OR WARRANTIES, EITHER EXPRESSED OR IMPLIED, OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR OF ANY OTHER NATURE ARE MADE HEREUNDER WITH RESPECT TO INFORMATION OR THE PRODUCT TO WHICH INFORMATION REFERS.

=====