

<p>SAREMCO AG Rohnacker CH-9445 Rebstein</p>	<p>EU Material Safety Data Sheet Microbond dentine conditioner</p>	<p>93/112/EWG page 1 of 4</p>
--	---	-----------------------------------

1. Commercial product name and supplier

commercial product name: microbond dentine conditioner

Application/use: conditioner

Manufacturer

SAREMCO AG
Rohnacker
CH-9445 Rebstein
Tel: 071 / 775 8060
Fax: 071 / 777 31 92

Notruf-Telefon: CH 071 / 775 8060

2. Composition

Chemical Characterization: Mixture of Acrylates, stabilizers, and initiators

Hazardous components:

- CAS Nr. 67-64-1 Acetone
F: highly flammable
- CAS Nr. 868-77-9 HEMA
Xi; R 36/38; R 43
- CAS Nr. 1830-78-0 Dimethacrylate
- CAS Nr. 110-16-7 Maleic acid
X_N; R 22; 36/37/38

3. Hazardous identification

Direct contact with uncured material can cause eye and skin irritation. The material is contraindicated if a person is known to be allergic to any of the ingredients of the product.

4. First Aid measures

- | | |
|---------------------|--|
| Inhalation | remove to fresh air |
| Eye contact | flush with plenty of water. Consult a physician if irritation persists |
| Skin contact | Wash thoroughly with soap and water |
| Ingestion | Consult a physician |
| Further information | if you feel unwell, seek medical advice; show this safety data sheet) |

5. Fire-fighting measures

Suitable extinguishing media: carbon dioxide, foam, dry chemicals

Extinguishing media to avoid: none known

Druckdatum: 25.5.2007	freigegeben: MH	Datum: 20.03.07	Ausgabe:07.01
-----------------------	-----------------	-----------------	---------------

<p>SAREMCO AG Rohnacker CH-9445 Rebstein</p>	<p>EU Material Safety Data Sheet Microbond dentine conditioner</p>	<p>93/112/EWG page 2 of 4</p>
--	---	-----------------------------------

6. Accidental release measures

Clean up mechanically; dispose of according to local and national regulations

7. Handling and storage

Handling: Only adequately trained personnel should handle this product. Keep out of reach of children

Storage: store in a cool and dry place. Avoid exposure to light

Fire and explosion protection: not required

8. Exposure controls/Personal protection

Technical measures: good ventilation should be sufficient

Control of treshold limits: none established

Personal protective equipment

Respiratory protection: not required

Hand protection; gloves

Eye protection: safety goggles

9. Physical and chemical properties

Appearance: liquid

Colour: yellowish

Odour: practically odorless

Chang of physical state -----

Melting point/melting area -----

Boiling point/boiling area: -----

Flash point: not determined

Ignition temperature: not determined

Lower explosion limits: not determined

Upper explosion limits: not determined

Vapour pressure: not determined

<p>Druckdatum: 25.5.2007</p>	<p>freigegeben: MH</p>	<p>Datum: 20.03.07</p>	<p>Ausgabe:07.01</p>
------------------------------	------------------------	------------------------	----------------------

<p>SAREMCO AG Rohnacker CH-9445 Rebstein</p>	<p>EU Material Safety Data Sheet Microbond dentine conditioner</p>	<p>93/112/EWG page 3 of 4</p>
--	---	-----------------------------------

Density not determined 20°C
Water solubility: < 0.1 %
pH-Value: < 7
Viscosity: not determined
Further information: none

10. Stability and reactivity

Thermal decomposition: None if used in accordance to instructions

Hazardous reactions: none

Hazardous decomposition products: none under normal conditions of storage and use

11. Toxicological information

Acute Toxicity: Oral LD50/ Ratte: 2000 mg/kg

Chronic toxicity: uncured material: prolonged or frequently repeated skin contact may cause allergic skin reactions

12. Ecologic information

No ecological problems to be anticipated if properly handled and used; nearly insoluble

13. Disposal considerations

The product should be disposed according to local and national regulations.

14. Transport information

Product is not classified for any mode of transportation

UN Number: ----

ADR/RID: ---

IMDG: ---

IATA: ---

<p>SAREMCO AG Rohnacker CH-9445 Rebstein</p>	<p>EU Material Safety Data Sheet Microbond dentine conditioner</p>	<p>93/112/EWG page 4 of 4</p>
--	--	-----------------------------------

15. Regulatory information

The product is a medical device according the EC-directive 93/42 EEC

EG-Nr.: ---

hazardous symbols: F highly flammable

risk phrases: R 11: highly flammable
R 22: harmful if swallowed
R 36/38: irritating to eyes and skin

Safety phrases: S 2: keep out of the reach of children
S 16: keep away from sources of ignition
S 26: in case of contact with eyes rinse immediatly with plenty of water and seek medical advise
S 28: after contact with skin wash immediatly with plenty of water and soap
S 24/25: avoid contact with skin and eyes
S 37/39: wear suitable gloves and eye/face protection

16. Other information

The above mentioned data correspon dto our present state of knowledge and experiance. The safety data sheet serves as a description of the products in regard to necessary safety measures. The indications do not have the meaning of guarantees on properties.